

CATALOG

PART I Safety Summary.....	2
PART II Product Use And Composition.....	3
PART III Main Specification Standard.....	3
PART IV Component Detail.....	
PART V Installation And Usage	4
PART VI Repair And Maintenance.....	5
PART VII After Sale Service	5

PART I Safety Summary

Warning :

- The user is strictly prohibited from being transported by any form of transport vehicle in a wheelchair.
- It is strictly forbidden to use the brakes during the driving. If there is no accompanying person, you can hold the handrim and gradually brake.
- It is strictly forbidden for the user to stand on the pedals, otherwise the wheelchair will fall over.
- The user is strictly prohibited from modifying the wheelchair without the consent of the technicians. company will not be responsible for any accident.
- Please do not stick out your body during driving, otherwise it will affect the stability of driving.

Notice :

- ① Practice using the wheelchair with a companion before using the wheelchair.
- ② On a wheelchair without an anti-roll device. Please do not sit on the wheelchair and do the rearward movement •
- ③ When using a wheelchair, please do not tilt the wheelchair to one side and do not lift the wheelchair.
- ④ Please park the wheelchair on a flat place.
- ⑤ Do not place tires near heaters.
- ⑥ When opening the wheelchair, please do not hold the seat tube and press down •
- ⑦ Adjustable pedals should be adjusted to a height of 50mm above the ground.
- ⑧ Please use it within the maximum weight range of the frame
- ⑨ Wheelchairs with anti-roll mechanism must be properly installed with anti-roll mechanism before use

We reserve the right to change the technology and appearance of this product. If there is any change without prior notice, please forgive me!


PART II Product Use And Composition

- It is used as a means of transportation for the disabled, the sick and the elderly and infirm.
- This product is mainly composed of frame, cross bar, armrest, armrest board, pedal, handbrake, hand grip, front wheel, rear wheel, brake, seat cushion, back cushion, etc.

PART III Main Specification Standard

1. Wheel landing: All wheels except anti-tip wheels must land smoothly.
2. Static stability: longitudinal $\geq 10^\circ$, lateral $\geq 15^\circ$.
3. Slope performance: $\geq 8^\circ$.
4. Sliding offset: $\leq 350\text{mm}$.
5. Minimum radius of gyration: $\leq 0.85\text{m}$.
6. Minimum conversion width: $\leq 1.5\text{m}$.

PART IV Component Detail



A	Front wheel
C	Seat cushion
E	Hand grip
G	Back cushion
I	Hand rim

B	Pedal
D	Armrest
F	Hand brake
H	Rear wheel
J	brake

PART V Installation And Usage

Open method

Hold the armrests with both hands so that the left and right frames are slightly separated. Press lightly on both sides of the seat cushion with your palms until the seat cushion is flat, and the wheelchair is ready to use.



Note: When unfolding, do not pull the armrest hard to avoid damage to various parts, and when pressing down on the seat cushion, do not hold your fingers on the seat tube to avoid pinching your fingers.

Folding method

First, turn up the left and right footrests, and gently lift the middle ends of the seat cushion with both hands to fold it.




Get on and off wheelchair method

•Get on wheelchair

When getting into the car, brake first, then hold the armrest tightly, and let your companion slowly help you into the seat.

• Get off wheelchair

Be sure to brake first, and then reverse the above method.



A	Front wheel
C	Seat cushion
E	Hand grip
G	Back cushion
I	Hand rim

B	Pedal
D	Armrest
F	Hand brake
H	Rear wheel
J	brake

PART VI Repair And Maintenance

After the wheelchair has been used for a period of time, it is necessary to check whether the movement of each part and whether the bolts are loose. If there is looseness, it should be adjusted in time to avoid unsafe hidden dangers.

A movable part should be lubricated weekly to prevent inflexible use.

- After using the wheelchair, use a soft dry cloth to wipe the surface moisture and pollution to prevent rust.
- Cushions and backrests should always be kept clean to prevent skin infection.
- The wheelchair should be placed in a dry, ventilated, stable and non-corrosive environment.
- Dispose of the waste in accordance with the city's regulations on environmental protection.

PART VII After Sale Service

This product has quality problems other than human factors within one week from the date of sale, and our company is responsible for refund, replacement and repair. Under normal use and storage, if the product has quality problems within one and a half years from the date of delivery (one year of use), the company will provide free maintenance. This product has quality problems after one and a half years from the date of delivery (one year of use). Users can go to the company's after-sales service department, office or dealer according to the invoice and warranty card, and the company will provide spare parts for maintenance. Reasonable charges. If the user is unable to provide the invoice, the warranty period shall be confirmed by the company's device number or ex-factory date extended by one month.

The following situations are not covered by the warranty:

- ① The front wheel and tire are worn or punctured, and the seat back cushion and armrest are damaged;
- ② Failures caused by unauthorized disassembly, repair or modification of the product;
- ③ Failure caused by accidental falling during use and handling;
- ④ Improper use, resulting in damage;
- ⑤ Failure to operate in the correct way in accordance with the instructions
- ⑥ Damage caused by unforeseen natural disasters.

