Elleboogkruk / Elbow Crutch / Ellenbogenkrücke

KL5804
KL5852**VEILIGHEIDSISTRUCTIES:**

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**SAFETY INSTRUCTIONS:**

Before using the crutches, it is advisable to have a specialist, or GP, look at them and help adjust them to your height. **CAUTION: The maximum user weight is 285 LBS (130 KG).** This must absolutely not be exceeded. Are the crutches already damaged before use? Then do not use them! Do not use the elbow crutches on wet and slippery surfaces. Do not hang bags or other objects on the crutches. This will negatively affect stability during use. Try to take into account the use of the crutches on your route as much as possible. Be extra careful on steep sections, stairs, thresholds and pavements.

SETTING THE HEIGHT:

Press the spring buttons of the leg and move up or down until the desired height is set. Then drop the spring knob into the hole until you hear a click. The handles should be at the height of the user's wrist when the user is in a natural, relaxed position with the arms along the body. Adjust the height of the cuff until it comes just below the elbow, and movement of the arm is not restricted.

CLEANING:

Clean the elbow crutches on a regular basis with warm soapy water. Do not use bleach. Always check that the chosen cleaning agent is suitable for cleaning the crutches. Do not use scouring pads or cleaning cloths with a rough surface as these products can permanently damage the crutches permanently, voiding the guarantee. Dry the crutches thoroughly before storing or using them again.

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