

KL5852

KL5804

EN

General guidelines

- Take the time to thoroughly read the following instruction and utilize the product accordingly.
- Ensure you retain this manual and pass it along when transferring the product.
- Please be aware that this summary may not encompass every detail of all variations and deliberated steps. Feel free to reach out to us for additional information and assistance.

Notes

- The maximum static load capacity is 130 kg (285 lbs), this must not be exceeded.
- Before assembly, carefully check the parts for damage. Improper assembly or damaged parts pose a risk of injury.
- Before using the crutches, it is advisable to have a specialist, or GP, look at them and help adjust them to your height.
- Under no circumstances, use the product if any part is missing or defective.
- Note: Do not use the elbow crutches on wet and slippery surfaces. You then risk slipping and injury.
- Do not hang bags or other objects on the crutches while using the crutches. This will adversely affect stability during use.
- Regularly check the strength of the joints.
- Try as much as possible to accommodate the use of the crutches on your route.
- Never use the elbow crutches for any purpose other than support.
- Take extra care when using the elbow crutches on steep sections, stairs, thresholds and pavements.
- This product is intended for indoor and outdoor use by one person.
- To set the correct height, the handles of the elbow crutches should be at the level of your wrist in a natural, relaxed position with the arms along the body.

Warnings

- Keep all parts and packaging materials (such as foil, plastic bags and foam) away from babies and children to avoid choking hazards.
- Keep small parts out of the reach of children during assembly. Swallowing or inhaling these parts can be fatal.
- Always use the elbow crutches correctly. It is strictly forbidden to use the elbow crutches as climbing aids.
- Children must be supervised by adults when using the product.
- Changes or modifications to the product without authorisation may compromise safety.

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