---

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---

| DUN

---

## General guidelines

Take the time to thoroughly read the following instruction and utilize the product accordingly.

Ensure you retain this manual and pass it along when transferring the product.

Please be aware that this summary may not encompass every detail of all variations and deliberated steps. Feel free to reach out to us for additional information and assistance.

## Warnings

Always consult a medical specialist before wearing support stockings, especially if you have health problems such as diabetes, severe circulation problems, heart failure, or inflammation in the legs.

It is essential to choose the right size support stockings for best effectiveness and comfort. Stockings that are too tight can lead to circulation problems, while stockings that are too loose lose their therapeutic effect.

**Note:** Do not wear the support stockings continuously. It is recommended to take them off at night, unless a doctor advises otherwise.

If you experience pain, tingling, swelling, skin discolouration or other unusual symptoms while wearing these stockings, seek medical attention immediately.

It is recommended not to wear the stockings over open wounds or untreated skin problems.

## Wearing Instructions

To enjoy maximum benefits, correct use of your Dunimed support stockings is essential.

1. Turn the sock inside out.
2. Place your foot in the toe piece first, then slip your heel in.
3. Lift the sock up over your calf, making sure it is evenly spread across your leg, to just below the knee.
4. Do not pull the sock further up or fold over the top edge.

## Washing instructions

To ensure that your Dunimed support stockings last as long as possible and maintain their optimal function, it is crucial that you wash and dry them properly. Here are detailed instructions on how to care for your socks:

**Washing method:** You can wash your Dunimed support stockings either by hand or in the washing machine. If you choose hand wash, do so gently in lukewarm water. When machine washing, it is important to set the temperature to a maximum of 30°C. This is essential to preserve the quality of the materials and prevent the socks from shrinking or losing their shape.

**Avoid tumble dryers:** These stockings cannot withstand the high temperatures and intense movements of a tumble dryer. Using a dryer can damage the fibres, causing the socks to lose their elasticity and support.

**Drying:** The best way to dry your Dunimed support stockings is to let them air dry. Hang them in a place with sufficient air circulation, but avoid direct sunlight. Direct sunlight can fade the colours of the stockings and affect the material. Make sure the socks are completely dry before storing or using them again.

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