Underwrap
Untertape

GLAD-4001

General guidelines

Take the time to thoroughly read the following instruction and utilize the product accordingly.

Ensure you retain this manual and pass it along when transferring the product.

Please be aware that this summary may not encompass every detail of all variations and deliberated steps. Feel free to reach out to us for additional information and assistance.

Warnings

Keep all parts and packaging materials (such as foil, plastic bags and foam) away from babies and children to avoid choking hazards.

Consult your doctor before using this product, especially if you suffer from skin problems, circulatory disorders or nervous disorders.

It is strictly forbidden to use the product for anything other than its intended use. Children should be supervised by adults when using the product.

Do not use the product if you have any disability or are taking any medication that may interfere with the safe use of this product.

In case of pain, tingling, swelling, skin discolouration or other unusual symptoms during use, stop using the tape and consult a doctor.

Do not use the tape over open wounds or untreated skin conditions. This may hinder the healing process or cause infections.

Notes

Make sure the skin is clean, dry and clear before applying the tape. In addition, make sure you apply the tape at room temperature; this ensures the best adhesion.

It is essential to avoid excessive tension when applying underwrap; a mild stretch is usually sufficient.

Long-term use of the tape without medical advice may have undesirable effects. Consult a doctor before long-term use.

Take extra care when used by children and the elderly. They may be more sensitive to the tape.

Store this product at room temperature in a dry, ventilated place to prevent mould growth.

Do not apply lotions, creams or ointments directly under the tape unless recommended by a doctor, as this may damage the tape or cause skin irritation.

When using the tape for sports activities, apply the tape at least 30 minutes before sports activities.

Replace the tape regularly every few days, depending on use. In case of sports activities, it is recommended to replace the tape immediately.

Remove the underwrap carefully to avoid skin damage, especially on sensitive skin.

Avoid wrinkles or creases in the underwrap as these can cause irritation.

Avoid using heating pads or ice directly over the underwrap, as this may damage the skin or affect the effectiveness of the tape.

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