Kinesiology Tape
Kinesio Tape

Glad-3929

General guidelines

Take the time to thoroughly read the following instruction and utilize the product accordingly.

Ensure you retain this manual and pass it along when transferring the product.

Please be aware that this summary may not encompass every detail of all variations and deliberated steps. Feel free to reach out to us for additional information and assistance.

Warnings

Keep all parts and packaging materials (such as foil, plastic bags and foam) away from babies and children to avoid choking hazards.

Consult your doctor before using this product, especially if you suffer from skin problems, circulatory disorders or nervous disorders.

It is strictly forbidden to use the product for anything other than its intended use. Children should be supervised by adults when using the product.

Do not use the product if you have any disability or are taking any medication that may interfere with the safe use of this product.

In case of pain, tingling, swelling, skin discolouration or other unusual symptoms during use, stop using the tape and consult a doctor.

Do not use the tape over open wounds or untreated skin conditions. This may hinder the healing process or cause infections.

Pregnant women should be careful when applying kinesio tape, especially on the abdomen and back.

Notes

Make sure the skin is clean, dry and clear before applying the tape. In addition, make sure you apply the tape at room temperature; this ensures the best adhesion.

It is essential to avoid excessive tension when applying kinesio tape; a mild stretch is usually sufficient.

Remove the tape carefully to avoid skin damage. It may help to remove the tape while showering.

Long-term use of the tape without medical advice may have unwanted effects. Consult a doctor before long-term use.

Take extra care when used by children and the elderly. They may be more sensitive to the tape.

Store this product at room temperature in a dry, ventilated place to prevent mould growth.

Do not apply lotions, creams or ointments directly under the tape unless recommended by a doctor, as this may damage the tape or cause skin irritation.

When using the tape for sports activities, apply the tape at least 30 minutes before sports activities.

Combining different types or colours of tape may be useful for specific therapeutic purposes.

Replace the tape regularly every few days, depending on use. In the case of sports activities, it is recommended to change the tape immediately.

When removing the tape, avoid pulling on the skin, which may cause discomfort or blisters.

Cut the kinesio tape to size with scissors, here rounding the corners is essential to prevent premature loosening of the tape.

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