



-1&1/F, Property 1, Area B, No.10, Lelin Road, Tongyi Industrial Park, Guzhen, Zhongshan, Guangdong, China

Phone:+86 0760-22323445

The lab is accredited for compliance with ISO/IEC17025:2017.

The results of tests, calibrations and/or measurements included in this document are traceable to Chinese/national standards.

CNAS is a signatory to the ILAC mutual recognition arrangement for the mutual recognition of the equivalence of testing, calibration and inspection reports.

TEST REPORT

Particulate respirator-half facepiece

EN 149:2001+A1:2009 Respiratory protective devices—Filtering half masks to protect against particles—requirements, testing, marking

The following samples were submitted and identified on behalf of the client as:

Product : Filtering Half Mask

Report No. : TSGK-2020-1887-T

Client : Universal Certification and Surveillance Service

Trade Ltd.Co.

Model(s) : XJ-9501

Date(s) of tests : 2020-09-27~2020-09-29

DESCRIPTION OF SAMPLES

General information

Classification
FFP2
White half mask

Manufacturer
FoShan Xujian Medical Equipment Co., Ltd

No.10 West Changxing Load, Junan district, Foshan City, Guangdong, China

A A A A A A A A A A A A

Approve: Reviewer: Chief Tester:

Issued:

陈健强

Test Report No. TSGK-2020-1887-T

Note:

The test results presented in this report relate to the samples tested only.

This report may be reproduced and distributed to your clients, provided that it is reproduced and distributed in full.

The authenticity of this test report and its contents can be verified by contacting thelaboratory.

Conclusion

Test Items

Clause 7.9.2	Penetration of filter material	Pass				
Clause 7.16	Breathing resistance	Pass				

Remarks: Pass = Meet EN 149:2001+A1:2009 FFP2 Requirement

Fail = Below EN 149:2001+A1:2009 FFP2 Requirement

N/A = Not Applicable

Disclaimer Measurement Uncertainty:

Unless otherwise agreed upon, Pass or Fail verdicts are given based on the measured values without any considerations of measurement uncertainties. Please note, every test method has a measurement uncertainty which has been evaluated by the laboratory according to ISO/IEC 17025 requirements.

By taking measurement uncertainties into account it might happen that measured values can neither be assessed as Pass nor as Fail.

Test Results

7.9.2 Penetration of filter material

 $^{\mathbf{l}}$

The penetration of the filter of the particle filtering half mask shall meet the requirements show in the table below.

Sodium chloride test 95 L/min

Paraffin oil test 95 L/min

FFP1	€20%	€20%
FFP2	≪6%	≪6%
FFP3	≤1%	≤1%

Note 1: FFP2 respirator. Test results are shown in Annex A Table 7.9.2.

7.16 Breathing resistance

Pass²

41	41	Maximum permitted resistance (mba	ur)
Classification	,cGF	Inhalation	Exhalation
	30 L/min	95 L/min	160 L/min
FFP1	0.6	2.1	3.0
FFP2	0.7	2.4	3.0
FFP3	1.0	3.0	3.0

Note 2: FFP2 respirator. Test results are shown in Annex A Table 7.16.

End of Test Results

Test Report No. TSGK-2020-1887-T

Annex A: Summarization of Test Data

Table 7.9.2 Penetration of filter material

Test specification: EN 149:2001+A1:2009 Clause 8.11

Aerosol	Condition	Sample No.	Penetration (%)	Assessment
	7	1	0.1	
Sodium chloride test	As received	2	0.2	
		3	0.2	
		4	0.6	Pass
Paraffin oil test	As received	5	0.8	
	7,	6	0.5	
Flow conditionin	ng: single filter: 95.0 L/min			

Table 7.16 Breathing resistance (mbar)

Test specification: EN 149:2001+A1:2009 Clause 8.9

	Flow rate		1			2				3							
As received			A	В	С	D	Е	A	В	C	D	Е	A	В	С	D	Е
	Inhalation -	30 L/min	0.5	0.5	0.4	0.4	0.5	0.4	0.4	0.4	0.5	0.4	0.4	0.4	0.5	0.4	0.5
		95 L/min	1.7	1.7	1.6	1.7	1.5	1.7	1.6	1.6	1.6	1.6	1.7	1.5	1.6	1.7	1.7
	Exhalation	160 L/min	1.2	1.2	1.2	1.1	1.1	1.2	1.2	1.1	1.2	1.2	1.1	1.2	1.1	1.1	1.2
Assessment	Pass																

A: facing directly ahead; B: facing vertically upwards; C: facing vertically downwards; D: lying on the left side; E: lying on the right side

End of Annex A

Annex B: photos of samples



End of Annex B